



Technology for Health!

Turkey's pride, the world's choice.
Novacheck, Blood Glucose Monitoring System.

We are thrilled to announce that we have been honored with the prestigious R&D Entrepreneur Award at the 13th biennial award ceremony organized by GAGIAD, one of Turkey's most influential and respected industrial associations.

This recognition is a testament to our commitment to innovation and excellence in research and development.

We are proud to share this significant achievement with our community and stakeholders!



MISSION

Our mission is to contribute to the protection and improvement of people's health by providing accurate and reliable diagnostic services to patients and healthcare professionals. To this end, we focus on reducing our environmental impact while working to increase the accessibility and quality of diagnostic services using the latest medical technologies.

This mission reflects the following core values: Adoption of the importance of diagnostic services for the protection and improvement of health; commitment to quality and ethical standards for the provision of accurate and reliable diagnostic services; respectful and dedicated approach to patients and healthcare professionals; and environmental sustainability and social responsibility awareness.

In line with this mission, our company will work to achieve the following goals:

- Provide the best diagnostic services to patients and healthcare professionals using the latest medical technologies.
- Work to increase the accessibility and quality of diagnostic services.
- Contribute to research and development in the field of diagnostics.
- Work to reduce our environmental impact.

Through these goals, our company believes that it will make significant contributions to the medical world and society.

Environmental Sustainability Initiatives: Efforts to reduce the environmental impact of facilities where diagnostic services are provided, reduce and recycle medical waste, and efficiently use energy and resources for diagnostic services reflect our commitment to environmental sustainability and social responsibility.

VISION

Novatech Diagnostics is a medical technology company dedicated to improving human health. We use our diagnostic products to help people live more active and fulfilling lives, feel better, live longer, and improve the quality of human life.

Our Vision:

To be a leader in diagnostic technologies to improve people's lives around the world. To continuously develop and innovate our diagnostic products and services. To be an innovative and pioneering company in the diagnostic industry.

Our Values:

Respect for Humanity: We view people as our most valuable asset and treat them with respect at all times.

Patient Focus: We always put patients' needs and wants first.

Transparency:

We always act in a transparent and honest manner in our business.

Integrity:

We always act in an honest and fair manner in our business.

Guided by this vision and values, Novatech Diagnostics is committed to improving human health and contributing to a healthier world.


BLOOD GLUCOSE MONITORING SYSTEMS NOVACHECK GLUCOSE METER

- Simple and intuitive design for ease of use
- Fast and accurate results in just 5 seconds
- Sensitive test strips that absorb blood quickly
- Large, bright display for easy viewing
- No coding required for quick and easy setup
- Tiny blood sample required for less discomfort
- Operates accurately at temperatures of 10-40° (50-104°)
- Operates in humidities of 10-85%
- Wide measurement range of 20-600 mg/dL (1.1-33.3 mmol/L)
- Stores up to 320 measurements for easy tracking
- Accurate measurement at hematocrit levels of 20-70%



NOVACHECK BLOOD GLUCOSE TEST STRIPS

- Glucose Dehydrogenase (GDH-FAD) enzyme for improved measurement accuracy
- Silver coating for effective and reliable results
- Capillary action for rapid and easy blood sampling
- Ergonomic design for easy insertion into the meter
- Durable construction to resist tearing and breakage



NOVACHECK PREGNANCY TEST

• Early Diagnosis, Quick Results, Reliable Performance

Novacheck pregnancy test offers the most advanced and user-friendly experience. With its cutting-edge technology and precise design, Novacheck pregnancy test assists in confirming pregnancy at early stages.

Advantages of Novacheck Pregnancy Test:

- High Sensitivity for Early Diagnosis: Thanks to advanced technology capable of detecting even low levels of hCG, it enables the early confirmation of pregnancy.
- Quick and Easy to Use: With its ergonomic design and simple instructions, the test can be easily administered in the comfort of your home.
- Results in Minutes: No need to worry for days waiting for results. Novacheck pregnancy test provides a clear answer within minutes.
- Reliable and Consistent Performance: High accuracy and repeatable results offer reliable information for healthcare professionals.
- Contributing to Early Diagnosis of Pregnancy: Novacheck pregnancy test helps protect the health of expectant mothers and their babies by facilitating early diagnosis.



%99.9
Accuracy
Doğruluk



%99.9
Accuracy
Doğruluk

- ◆ Fast
- ◆ Easy
- ◆ Reliable
- ◆ Hygienic



NOVACHECK SARS CoV-2 RAPID ANTIGEN TEST

Novatech Diagnostics is the exclusive manufacturer in Turkey listed in Category A of the European Commission's Common Rapid Antigen Test (RAT) List.



Prospective study
from **Heidelberg**, Germany

Quality
tested in
Germany

One of the best field study results
Novacheck[®] Antigen Rapid Test




Omicron detection
100% of BA.1 and BA.2

One of the best field study results Antigen Rapid Test



Novacheck®

Prospective study
from Heidelberg, Germany



Ct-Value	Number of the Samples	Number of true positive AG rapid test samples	Number of false negative AG rapid test samples	Sensitivity of the AG Rapid Test
≤ 25	60	59	1	98.36%
≤ 30	82	81	1	98.78%
≤ 32	92	87	5	94.57%
≤ 34	101	91	10	90.10%
≤ 37	106	92	14	86.79%
Overall*	110	93	17	84.55%

*including samples with a Ct value > 37

Method	Real-time PCR			
	Result	Positive	Negative	Final result
Novacheck	Positive	91	6	97
	Negative	10	357	367
Final result		101	363	464


Sensitivity: 98.36%
Specificity: 98.30%
Accuracy: 96.55%
Ct≤25



EU Common RAT list of COVID-19
Antigen tests A-Category



Omicron detection 100% of BA.1 and BA.2





Ct-Value	Number of the Samples	Number of true positive AG rapid test samples	Number of false negative AG rapid test samples	Sensitivity of the AG Rapid Test
≤ 30	7	7	0	100%
≤ 32	18	18	0	100%
≤ 34	30	30	0	100%
≤ 37	41	41	0	100%


Diagnostic sensitivity with respect to the Omicron variant of SARS-CoV-2 100%






 www.novadiag.com

 info@novadiag.com

 +90 342 502 27 27

 2. Organize Sanayi Bölgesi Hacı Sani
Konukoğlu Bulvarı 83228 Nolu Cad. No:17
27620 Şehitkamil / Gaziantep / Türkiye



website



Technology for Health!